

# MC Disinfectant Wipe

## Pre-moistened Wipes

### TECHNICAL BULLETIN

#### Where to Use

This product is intended for use in:

Airports, Funeral homes, Barber shops, Daycare centers, Gyms, Athletic facilities, Beauty salons, Hospices, Hospitals and healthcare settings: Ambulances, Anesthesia, CAT labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's office(s), Emergency Rooms, E.R., Home health care, ICU, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating rooms, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physician's office(s), Ophthalmic offices, Radiology, Recovery Room, Rehabilitation, Respiratory therapy, Transport vehicles.

#### Surface Safe

This product will clean and disinfect hard, non-porous, surfaces of: bed railings, cabinets, carts, chairs, counters, exam tables, gurneys, isolettes, IV poles, stethoscopes, stretchers, tables, telephones, toilet seats, toilets, bathtubs, bathroom fixtures, shower stalls, sinks, tiles, vanity tops, hospital beds, traction devices, MRI, CAT, scales, paddles, wheelchairs, ultrasound transducers and probes.

This product will clean and disinfect hard, non-porous surfaces of: grocery cart handles, grocery carts, shopping carts, grocery cart child seats.



Item Number 953398PL

### Multi-Clean MC Disinfectant Wipes

(DIN No. 02523280)

Kills SARS-CoV-2, the cause of COVID-19  
on hard, non-porous surfaces in just 30 seconds!

#### ACTIVE INGREDIENTS

Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides .....	0.14%
Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides .....	0.14%
INERT INGREDIENTS .....	99.72%
Total .....	100.00%

DIN No: 02523280

MC Disinfecting Wipes are intended For Use in Domestic, Institutional/Industrial, Hospital, Food Processing and Barn Premises: ambulances, anesthesia, CAT labs, CCU, Central Supply, Clinics, Dental Offices, Dialysis Clinics, Doctor's office(s), Emergency Rooms, E.R., Home health care, Hospices, ICU, Laboratory, Laundry rooms, Neonatal Intensive Care Units (NICU), Newborn nursery, Nursing homes, Operating rooms, Pediatrics, Pediatric Intensive Care Units (PICU), Physical therapy, Physician's office(s), Ophthalmic offices, Radiology, Recovery Room, Rehabilitation, Respiratory therapy, Surgery, Transport vehicles, X-Ray, Animal holding areas, veterinary waiting rooms, veterinary examination rooms, veterinary X-ray rooms, veterinary operating rooms, animal transportation vehicles.

**DIRECTIONS FOR USE:** This product is not for use on critical and semi-critical medical device surfaces.

**CLEANING:** Use one MC Disinfecting Wipes to completely remove all loose debris from surfaces to be treated.

**DISINFECTING:** To disinfect hard, non-porous surfaces, thoroughly wet the surface to be treated. Treated surface must remain visibly wet for disinfecting. Use enough wipes for treated surface to remain visibly wet for 5 minutes (30 seconds for SARS-CoV-2). For heavily soiled surfaces, clean before following disinfecting instructions. Avoid contact with food. Remove all animals/poultry and their feed from premises prior to disinfection.

Remove all heavy soil, such as urine and fecal matter, from surfaces and objects prior to disinfection. Empty all feeding and watering appliances prior to disinfection. Following disinfection, do not house animals/poultry until areas have been ventilated. All surfaces and/or objects that will contact feed or drinking water should be rinsed with potable water before reuse.

MC Disinfecting Wipes is an effective disinfectant against *Escherichia coli* (E. coli) O157:H7, *Pseudomonas aeruginosa*, *Salmonella enterica* (Salmonella), *Staphylococcus aureus*, *Streptococcus pyogenes* (Strep), Community Associated Methicillin Resistant *Staphylococcus aureus* CA-MRSA (Genotype USA400) and Community Associated Methicillin resistant *Staphylococcus aureus* CA-MRSA (Genotype USA300), Vancomycin resistant *Enterococcus faecium* (VRE), Methicillin resistant *Staphylococcus aureus* (MRSA), Vancomycin intermediate resistant *Staphylococcus aureus* (VISA), *Klebsiella pneumoniae* New Delhi Metallo-Beta Lactamase (NDM-1) Carbapenem Resistant.

MC Disinfecting Wipes is an effective disinfectant on hard, non-porous surfaces against Influenza A2/Japan, Herpes Simplex Virus 1 (ATCC VR-260) and Herpes Simplex Virus 2 (ATCC VR-734) when the treated surface is allowed to remain wet for 5 minutes (30 seconds for SARS-CoV-2 the virus that causes COVID-19).

This product has demonstrated effectiveness against Influenza A virus and is expected to inactivate all Influenza A viruses including 2009 (H1N1) pandemic Influenza A virus.

**FUNGICIDAL ACTIVITY:** MC Disinfecting Wipes is an effective fungicide against *Trichophyton interdigitale* (ATCC 9533), when used as directed on hard surfaces found in bathrooms, shower stalls, locker rooms, exercise facilities or other clean, hard, non-porous, surfaces commonly contacted by bare feet. The surface must remain visibly wet for 3 minutes. Allow surfaces to air dry.

**EFFICACY TESTS HAVE DEMONSTRATED THAT THIS PRODUCT IS AN EFFECTIVE BACTERICIDE AND FUNGICIDE IN THE PRESENCE OF ORGANIC SOIL (5% BLOOD SERUM).**

**SANITIZING/DEODORIZING:** To sanitize /deodorize hard, non-porous surfaces, wipe surface, using enough wipes to allow treated surface to remain wet for 60 seconds for *Staphylococcus aureus* (Staph), *Klebsiella pneumoniae*, *Listeria monocytogenes*, *Escherichia coli* (E. coli) O157:H7, *Campylobacter jejuni*. Let surface dry. For highly soiled surfaces, clean excess dirt first.

For surfaces and/or objects that may come into direct contact with children at the mouthing stage of development, a rinse with potable water is recommended.

For surfaces that may come in contact with food, a potable water rinse is required.

For plastic and painted surfaces, spot test on an inconspicuous area before use.

MC Disinfecting Wipes will not harm most surfaces including acrylic, fiberglass and vinyl.

Not recommended for use on unpainted wood.

Not for personal cleansing.

This is not a baby wipe!

Do not flush in toilet.

**PRECAUTIONARY STATEMENTS HAZARDS TO HUMANS AND DOMESTIC ANIMALS CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling.

**FIRST AID:** Have the product container or label and Drug Identification Number with you when calling a poison control center or doctor, or going for treatment.

**IF IN EYES:** Hold eyes open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

## MICRO-ORGANISMS

### Viruses:

Bovine Viral Diarrhea Virus (BVDV).....	5 minutes
Influenza A2/Japan/305/57 Virus (H2N2).....	5 minutes
Hepatitis B Virus (HBV).....	5 minutes
Hepatitis C Virus (HCV).....	5 minutes
Duck Hepatitis B Virus (DHBV).....	5 minutes
Herpes Simplex Virus 1 (HSV-1).....	5 minutes
Herpes Simplex Virus 2 (HSV-2).....	5 minutes
Human Coronavirus.....	10 minutes
Human Immunodeficiency Virus Type 1 (HIV-1).....	5 minutes
SARS-Related Coronavirus 2 [SARS-CoV-2]	
USA-WA1/2020 BEI NR-52281.....	30 seconds

---

### Bacteria:

Pseudomonas aeruginosa.....	5 minutes
Staphylococcus aureus.....	5 minutes
Salmonella enterica (Salmonella).....	5 minutes
Vancomycin Resistant Enterococcus faecium (VRE).....	5 minutes
Escherichia coli O157:H7.....	5 minutes
Methicillin resistant Staphylococcus aureus (MRSA).....	5 minutes
Streptococcus pyogenes.....	5 minutes
Vancomycin intermediate resistant Staphylococcus aureus (VISA)...	5 minutes
Klebsiella pneumoniae New Delhi Metallo-Beta Lactamase (NDM-1) Carbapenem Resistant.....	5 minutes
Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA) Genotype USA300.....	5 minutes
Community Associated Methicillin resistant Staphylococcus aureus (CA-MRSA) Genotype USA400.....	5 minutes

---

### Fungicidal:

Trychophyton interdigitale, formerly known as trichophyton mentagrophytes	3 minutes
--	-----------



600 Cardigan Road  
Shoreview, Minnesota  
www.multi-clean.com